



PHYSICS

UNIVERSITY OF TORONTO

Book #2

2025-02

Subject ADVANCED PHYSICS LAB

Name ADITYA K. RAO

Address ADI.RAO@MAIL.UTORONTO.CA

ADITYA KANURU RAO
1008307761
ADI.RAO@MAIL.UTORONTO.CA
+1(416) 802-5648

SUMMARY ABSTRACT

Nuclear Magnetic Resonance is a spectroscopic technique used to determine the structure, composition, and dynamics of nuclei. It is based on the principle that the magnetic constituents of such nuclei are driven by a magnetic moment and an angular momentum component (the latter of which is quantized). In this experiment a TechSpin NMR Spectroscope to measure various properties of Heavy Mineral Oil, Light Mineral Oil, and Distilled Water.

Namely, investigations focused on tuning the machine, measuring the T_1 , spin-lattice relaxation time, the FID signal, and attempting to measure the spin-spin relaxation time known as T_2 . Additional exploration into the applications of NMR in other fields was also conducted, namely, Quantum Computing where NMR is used as techniques can be used to implement single and two qubit unitary gates in Trapped Ion systems.

Experimental findings yielded $T_1 \approx 447.8 \pm 0.8 \mu\text{s}$, though this may be flawed as a proper resonance condition was not achieved. Note that this measurement was for mineral oil (heavy).

T_2 and T_1 (cyclic) measurements were made, but error in data collection meant ~~adequate~~ insufficient data was obtained to properly calculate.

Additions to the experiment include a Foray into Quantum Computing, Measurements on Spin-Spin diffusion (T_2) for H_2O (distilled) and a non-inverting operational amplifier circuit being made to boost low signals.

Future endeavors would primarily be guided by better time management and more complete data collection.

Acknowledgements :

I would like to thank Professor Richard and Teaching Assistant Kyle Thompson for their guidance when conducting the lab. Additionally, I would like to thank Professor Aria Harik and Professor Tahir Shaban for their ideas on extensions for the lab. Finally, I would like to thank the lab technicians for helping me remove the distilled H_2O sample when it almost got stuck in the machine.

TABLE OF CONTENTS

Page Range	Title of Section
01	Summary Abstract
02	Table of Contents
03	Quick Reference Material
04	References
05 - 09	Background Theory
10 - 12	Experimental Planning
13 - 15	Equipment Manual Notes
16 - 19	Equipment Setup
20 - 20	Equipment Calibration
21 - 24	Initial Observations
25 - 28	Data Collection For Heavy and Light Mineral Oil
29 - 30	Data Analysis For Heavy and Light Mineral Oil
31 - 33	Data Collection For Spin Echo in Heavy and Light Mineral Oil
34 - 37	Adjustment of Resonance Finding Methodology
38 - 39	New Data Collection for Heavy and Light Mineral Oil
40	Beginning T ₂ Measurement in Distilled Water
41 - 42	Design of Non-Inverting Op-Amp Circuit for Measurements
43 - 47	Measurements for T ₂ in Distilled Water
48 - 50	Exploration into Unitary Gates
51 - 53	Literature Review: Application of NMR in Quantum Computing
54 - 55	Final Comments on Analysis + Mini Reflection on Background Experiments
55 - 71	Appendix (Code + Data)

QUICK REFERENCE MATERIAL

Procedures

Pg #	Code
14	PO2.N1A
16	PO2.N1B
18	PO2.N1G

Purpose

Experiment Apparatus Connection Setup Verification
 Experiment Apparatus Cable Connections
 Single Pulsed NMR Experiment Procedure

Key Contacts

~~NAME~~

Name: Prof. Pierre Savard
 Email: SAVARD@PHYSICS.UTORONTO.CA

} Supervising Professor

Name: TA. Kyle Thompson

Email: KTHOMPSON@PHYSICS.UTORONTO.CA

} Supervising Teaching Assistant

Name: Prof. Anja Harlick

Email: ANJA.HARLICK@UTORONTO.CA

} Course Coordinator
 (PHM426)

Personal Contact Info

← Should this book be found, please contact me at
 one of the following methods.

Name: Aditya Kanuru Rao (Adi.Rao)
 Email: ADIRAO@MAIL.UTORONTO.CA
 Phone: +1 (416) 802-5648
 GitHub: @ADIRAO-PROJECTS

ADITYAKRAO.PROJECTS@GMAIL.COM

REFERENCES

- 1 D. Bailey, J. Borwein, S. Vire, and J. Pitre, NMR Nuclear Magnetic Resonance, edited by J. Shepard (University of Toronto, 2013)
- 2 TeachSpin Inc., TeachSpin Inc. Pulsed Spin NMR Spectrometer Instructor's Manual
- 3 A. Trabesinger, "Long live the spin," Nature Physics 8, 781-782 (2012)
- 4 T. Wogan, "NMR Spectroscopy without the "M"- physics world" (2011)
- 5 L. Yanns, "Researchers increase NMR/MRI sensitivity through hyperpolarization of nuclei in diamond", (2012)
- 6 L. Zyga, "143 is the largest number yet to be factored by a Quantum Algorithm" (2012)
7. M.A. Nielsen and I.L. Chuang, Quantum Computing and Quantum Information (Cambridge University Press, 2010)

25.02.09 BACKGROUND THEORY
Nuclei

3.47 Nuclear Magnetic Resonance, or NMR, is a phenomena primarily observed in systems where the magnetism ^{couplings} has both a magnetic and angular ^{momentum} property. Almost all stable nuclei have this

51 classically, we may consider a spinning dipole to emulate such nuclei:




Fig 02.09 A.1

,59 which induces this magnetism. Then, we may link the angular momentum \vec{J} with the gyromagnetic ratio γ to obtain the resultant magnetic moment $\vec{\mu}$

unique for each nucleus

$$\vec{\mu} = \gamma \vec{J}$$

406 From Quantum Mechanics, we also know that the angular momentum is indeed quantized leading to a direct relation b/w angular momentum and the spin of the particle

$$\vec{J} = (\hbar S) \text{--- spin}$$

$\frac{\hbar}{2\pi} \leftarrow \text{Planck's constant}$

508 We may then calculate the total magnetic energy in an external field as

$$\vec{E}_M = -\vec{\mu} \cdot \vec{B}$$

522 In our experiment, we will primarily only be applying this in the z direction leading to a simplification of the equation to

$$E_M = -\mu_z B = -\gamma \hbar \vec{S} \cdot \vec{B}$$

Quantized w/

$$\vec{M} = \vec{S} \rightarrow S^2 = n$$


527 For fermionic spin- $\frac{1}{2}$ systems, this leads (\pm spins) this leads to the magnetic allowed values of S to be $\pm \frac{1}{2}$ or spin up/down in the up respectively. (I may use $|+\rangle$ and $|-\rangle$ to denote these in the unpolarized basis).

1532 Then our allowed energies become

$$E_{M_0} = -\gamma \hbar (-\frac{1}{2}) B_0 = \frac{\gamma \hbar}{2} B_0$$

$$\text{and } E_{M_1} = -\gamma (\pm) (\frac{1}{2}) B_0 = -\frac{\gamma \hbar}{2} B_0$$

1539 which may be visualized as follows



1611 which implies our energy can be written in terms of some angular frequency $\omega_0 = \gamma B_0$

$$\Rightarrow \Delta E_M = \hbar \omega_0$$

1613 Note that this only is in a constant magnetic field (B_0) in the \hat{z} direction.

From Eq. 2, we may write for a given:

$$\gamma_B = 2.675 \cdot 10^8 \text{ rad (ST)}^{-1}$$

1616 and we get $f_{\text{res}} = 42.58 B_0 \text{ T}$ resonant frequency.

Then, say we have N_0 spins in the downstate and N_1 in the up state in Thermal Equilibrium, then the population ratio of the spins is given by the Boltzmann factor

$$\frac{N_1}{N_0} = \exp \left\{ \frac{\Delta E_M}{kT} \right\} = \exp \left\{ \frac{\hbar \omega_0}{kT} \right\}$$

with magnetization

$$M_z = (N_0 - N_1) \mu$$

1627 However, in our NMR experiment, this will not appear instantly. Instead the Spin relaxation time (Spin-Lattice Relaxation time) known as T_1 is given from solving the ODE

$$\frac{d_t M_z}{dt} = \frac{M_0 - M_z}{T_1}$$

$$M_0 = (N_0 + N_1) \tanh\left(\frac{\mu B}{kT}\right)$$

$$\approx N \frac{\mu^2 B}{kT}$$


16.37 Placing an unmagnetized sample implies that at $t=0, M_z=0$,

$$M_z(t) = M_0 \left(1 - \exp\left\{-\frac{t}{T_1}\right\}\right)$$

which can be used to find our T_1 time later on.

16.43 Pulsed NMR (what this experiment deals with) is especially suited to measuring such relaxation times as it can make "direct unambiguous measurements" [6]

16.45 However, our machine may only measure in the xy plane. This can be understood classically



16.58 Then, through some manipulation, we may get the net magnetization in the xy plane as for non-constant magnetic field.

$$d_t \vec{M} = \gamma \vec{M} \times \vec{B} \quad \text{and} \quad \vec{B}(t) = \underbrace{B_0 \cos \omega t \hat{x} + B_0 \sin \omega t \hat{y}}_{+ B_0^2} + B_0^2$$

17.20 We may manipulate these fields to manipulate the magnetism. If we create a rotating field at our resonant frequency

$$\omega_0 = \gamma B_0$$

17.25 The result is that the magnetization begins to precess about the driving field at $\omega = \gamma B_1$, as can be seen below.




Fig 02.09c.2

17.35 Depending on the duration of this magnetic pulse, the following transformations will be observed.

$$\frac{\pi}{2} \longrightarrow M_z \longrightarrow M_y \quad (\text{Author Rotation})$$

$$\pi \longrightarrow M_z \longrightarrow -M_z \quad (\text{Half Rotation})$$

$$2\pi \longrightarrow M_z \longrightarrow M_z \quad (\text{Full Rotation})$$

17.43

2025-02-12

1609

~~Signaling out~~ Note this bit of background was written 2025-02-11 so I have omitted some graphs. Given that an aluminized spectrometer only measures in the M_y plane, we may see only some of the magnetization components as it passes in the M_y plane.

16.15 We may solve for this as

$$\partial_t M_{x'} = -\frac{M_{x'}}{T_2} \quad \text{AND} \quad \partial_t M_{y'} = -\frac{M_{y'}}{T_2}$$

w/ solutions given by

$$M_{x'} = M_0 \exp\{-t/T_2\} \quad M_{y'} = M_0 \exp\{t/T_2\}$$

where T_2 is known as the spin-spin relaxation time.

1626 which we may record by seeing how long it takes a spin to return to its equilibrium state. which may look as the following:




Fig 02.09 D.1 : sketch of FID (Free Induction Decay)
Signal expected for spin-spin relaxation w/ influence from magnetite.

~~lett~~

1901 In the figure, this may look like the T_2 relaxation time, but is actually the Free Induction Decay (FID) which is actually due to the non-uniformity of the magnetic field. For the PS2 controller [2] if the $T_2 < 3\text{ ms}$, then the FID $\propto T_2$ time since it is accurate enough to produce a magnetic field of uniformity of 0.3 ms decay time which can be improved using shims coils to improve the T_2^* time (decay time due to magnet).

1911 Even without these methods, we may use E. Harris's discovery of a spin-echo to still accurately measure T_2 [2]. By comparing the maximum height of the Echo and FID signal, one would notice a distinct difference.

1932 This loss occurs due to "fluctuations at the nuclear sites" [2] which are a measure of T_2 . Then, varying the delay b/w signals, T_2 may be obtained from the equation

$$M_{xy}(2\tau) = M_0 \exp \left\{ -\frac{2\tau}{T_2} \right\}$$

1947 Signing out ~~AFK~~

10

2025-02-04

1029

Signing into lab

1030

Setting up lab book

1101

Discussions w/ Prof PS and TA KT

- Make sure to use o-rings before putting sample into magnet
- Document figures well

1200

Figure setups done, reading lab manual & literature.
+ setting up filestructure and code base

1320

Setup preliminaries & examination setup.

1358

Not much completed today other than cleaning station and getting set up.
Likely will review manual & background today and ~~attempt~~ and
plan out Friday session (+ roadmap for full ~~one~~ experiment).

1400

Signing out of lab *Auf Wiedersehen*

11/20/07

11/21

11/22

11/23

Signing into lab

Completing some grading work for 11:00, did some light review for NMR last night.

Schematic of Spectrometer




Figure 02.07 A.1 : Schematic of Spectrometer

11/24 As in Fig 02.07 A.1, we place the sample w/ the pickup coils

11/25 Inspecting apparatus

11/26 Setup below





FIGURE 02.07 A.2: Initial experimental setup for NMR.



T202 Experiment Plan & Reading Literature

1. Initial Measurements
2. Fundamental Literature Review
3. Advanced Measurements
4. Data Analysis
5. Extension
6. Literature Review in Quantum computing or other field.
7. Literature Review in Frustrated magnetic States.

13.09 Signing out of lab AR

5-02-11

910 Signaling in NMR

0916 Objectives

1. Write in Procedures
2. Some data collection

0920 From the F02.07A.2, the sample sits w/i the coil. B/w the transmitter and receiver, notice a pair of crossed diodes. These restrict the voltage from the sample to ~0.6V (reduction of noise).

The crossed diodes to ground protect input from overvolt from transmitter
 \Rightarrow NMR output in mV

when far from Resonance, we expect to see Fig 02.11A.1 at the Probe Head and Fig 02.11A.2 at LF-HF Output w/ the RC Filter at HF.




FIGURE 02.11A.1: Conceptual sketch of r.f. pulses at the PROBE HEAD from the lab manual [1].




FIGURE 02.11A.2: r.f. pulses at the LF-HF Output for RC Filter set to HF from the lab manual [1].

0935

One experiment we will conduct as a preliminary is finding the resonance of FeCl_3 .

0936

Procedure (P02.11A)

1. Set RC filter to 0.5
2. Connect LF-PAF out to oscilloscope
3. Verify settings correspond w/ obs output.
4. Verify Magnet coolant is ON
5. Verify Magnet Current & Voltage controls set at zero
6. Set Magnet PSupply to voltage limiting mode { AS temp 1
Voltage & mag drift
 - (6.1) Set current to Max (Voltage ZERO)
 - (6.2) Turn up voltage until induction coil observed
DO NOT EXCEED 7A
 - (6.3) Turn off Reverse step 6
7. To Turn off magnet, reverse step 7 \rightarrow 1

0942 Performing Proc. P02.11A

0942

Ready Lanthanum Ion Tech Spintron

0954

Summary of key parts

Full Benchmarking will be evaluated later

$$\text{Resonance freq } \omega_0 = \gamma B_0$$

$$\gamma_{\text{polar}} = 2.675 \cdot 10^8 \text{ rad s}^{-1} \text{ T}^{-1}$$

 $\frac{d\theta}{dt}$

$$f_{\text{polar}} = 42.58 \text{ MHz T}^{-1}$$

$$f_{\text{Lanthan}} = 40 \cdot 0.055 \text{ MHz T}^{-1}$$

Thermal equilibrium // H field, approaches eq. behavior
exp rate w/ T_1 spin-弛豫 time.

Classically

$$\frac{d}{dt} \vec{\mu} = \gamma (\vec{m} \times \vec{B})$$

May be time
dependent

→
57 Pulsed NMR employs rotating RF described by

$$\vec{B}(\omega) = B_1 \cos \omega t \hat{x} + B_1 \sin \omega t \hat{y} + B_0 \hat{z}$$

$$\vec{H} = B_1 \hat{x}' + (B_0 - \omega/k) \hat{y}' \leftarrow \text{Effective field}$$

when on resonance $\omega = \omega_0 = \gamma B_0$ and $H = B_1 \hat{x}'$ and precesses around $\vec{B}_0 \cdot \vec{H}$.

$\pi/2$ Pulse : left on long enough for Mo to rotate about $\hat{x}\text{-}\hat{y}$ plane.

$$t_w = \frac{\pi}{2\omega_1} \quad \text{w/ } \omega_1 = \gamma B_1 \\ \Rightarrow t_w(\pi/2) = \frac{\pi}{2\gamma B_1}$$

T_2 is spin relaxation time

Spin echo allow measurement of T_2 shown below in Fig 02.11B.1




Figure 02.11 B.1: Spin Echo Generation using $\frac{\pi}{2}$ -pulse and π -pulse to measure spin relaxation time T_2 of sample.
Figure sketch from [2]

1011 Notes on Instrument from "Getting Started" in [2]



1014

Procedure P02.11B

1 Connect the Blue 18in cable from Synth Pulsed RF out to Pulsed RF In Rce

2 Connect the following BNC connectors

Q(PP) w/ Q(Synth)

I(PP) w/ I(Synth)

Blanking out (PP) w/ Blanking out (Synth)

RF out (Synth) w/ RF In (Rce)

3 Connect the following to scope.

Sync Out (PP) w/ Trigger

Env. Out (Rce) w/ CH1

Q/I out (Rce) w/ CH2

4 Perform / Verify the following

CW out (Synth), B Pulse(PP), M4(PR) OFE

Sync to A (PP)

Toggle

Pulse to A (PP)

Toggle

Ret out (Synth) ON

Filter TC to 0.1 (Rce)

Train 15%

Band to 8 (Rce)

Blanking on (Rce)

width 25% (Rce)

Do NOT CONNECT

Rce (Rce) CW IN RF OUT

Synth (Synth) RF out SWEEP IN, CW OUT

PulseProg (PP) EXT start

lock-in All CONNECTORS

1023

Verifying Proc. P02.11B

1023

Verifying Do Not Connect

CW IN ✓

SWEEP IN ✓

Ext Out ✓

RF OUT ✓

CW OUT ✓

lock In (all) ✓

→ Verifying ~~task~~ (P02.11B-1) ✓

Verifying ~~task~~ (P02.11B-2)

$$\begin{aligned} Q(\text{PP}) &= Q(\text{Synth}) \quad \checkmark \\ I(\text{PP}) &= I(\text{Synth}) \quad \checkmark \end{aligned}$$

$$\begin{aligned} \text{Blanker (PP)} &= (\text{Synth}) \quad \checkmark \\ \text{RF out (Synth)} &= \text{Ref In (RF)} \quad \checkmark \end{aligned}$$

Verifying (P02.11B-3)

(Sync out (PP) w/ trigger)

Digital Oscilloscope - DSOX1204C has no trigger input, will perform manual / auto triggers

1030

Env out w/ ch1 ✓

Q/I out (Rec) w/ ch2

Q/I out { Switched Q to ch2
on Receiver { Switched & I to ch3

1031 Verifying (P02.11B-3)

CW out, B Pulse, MGT
OFF ✓ OFF ✓ OFF ✓

~~Sync~~ Sync Pulse
A ✓ toggled to A ✓

Ret out
ON

Gain
75%
Blanking
ON

FILTER T/C
Adj to 0.1

Band
?

Width

{ No gradations
in width supply
set to 25%.

1032 No gradations on width nob, ~~sync~~ roughly read just to 75%. Used masking logic and sharpie to mark gradations.

1033 Set up completed his connections

1045

Discussion w/ Prof DS

1. Need to calibrate.
2. Needed to explain NMR and taught the physics at the experiment
3. Will go back at end of session.

→ 1046 Procedure (P02.NC)

Module Parameters adjusted w/ Path to set rising encoder.

Syst Default settings | Pulse Programmer Default

F Frequency 18.00000 MHz
 P Refr Phase -180°
 A CW Power -10 dBm
 S Sweep 0 kHz V⁻¹

A A length 0.02 μs
 B B length 0.02 μs
 T Tau 0.0001 s
 N Number B 0
 P period 0.2 ms

Lock-In Default

G Gain 80V
 P Reference Phase -180°
 T Time const. 0.5s
 M Modulation Amp - OFF
 H Field off 0.00G
 S Sweep Modul OFF
 A Sweep Amp. 0.0314
 D Sweep Duration 10s

1052 Procedure (P02.NC) : Single Pulse Experiment for Mineral oil

will take lower miniscus

- 1 Ensure sample only fills 3mm of tube
- 2 Load Oring Stop at 39mm
- 3 While Picking Probe to measure RF fields during pulse
Useage maximum at 65mm
- 4 While Picking Probe to tune to Larmor Frequency at bottom
- 5 Place probe into sample solenoid
- 6 Set Syst freq to 18.00
- 7 Set A length to 2.5μs
- 8 Set P to 100mW
- 9 Set Sweep function
- 10 Subtract Min off max

101 Discussion w/ TA KT

→ Turned on Machine

→ Explained how to tune capacitors of magnet

102 Inserted Pickup Probe into sample area, plugged probe into ch4. Adjusting capacitors with maximum ohm meter.

103 Turning off equipment. Noticing error signals on PS2 controller

106 Tuning capacitor Schematic Fig 02.11C-1




Figure 02.11C-1

109 Turning on equipment and adjusting to default points.

110 Verified Rebut settings. Changing rebut settings to those in PROC. (PO211C - 6-8)

Synth 20.90000 MHz

A-tens to 2.5 μs

P to 100 ms

114 Opting to change scope on the fly.

115 Pickup probe placed in sample.

117 Scope image in Fig 02.11 D.1




FIGURE 02.11D.1: Initial Reading from Oscilloscope. Notice Peaks

1118 Adjusting capacitors until 40Vpp observed.

1122 Only observing maximum V_{pp} of 32 V

1133 Discussion w/ TA, KT and LT, CA

$\sim 33 V_{pp}$ should be okay despite what the manual says.

Going to continue at next meeting. Images others




FIGURE 02.11D.2: Scope readings using zoom mode. Clearly we can see the peaks.




FIGURE 02.11D.3: Another scope reading using zoom mode. Maximum voltage observed $V_{pp} = 31.5 \pm 0.5$ V

1204 Preparing mineral oil samples for Free Induction Decay measurement

Seeing as this is just a test, will use preprepared sample.

1210 Connect CH1 to ENV out, disconnected & removed pickup probe

I and Q outputs are the product of the signal from precessing spins with the reference signal from the oscillator

$$\sin(\omega_{\text{rf}} t) \sin(\omega_{\text{spins}} t) = \frac{1}{2} \cos(\underbrace{\omega_{\text{rf}} t - \omega_{\text{spins}} t}_{\text{Beat signal}}) - \frac{1}{2} \cos(\underbrace{\omega_{\text{rf}} t + \omega_{\text{spins}} t}_{\text{Difference signal}})$$

Should be filtered out ~40MHz

Result at output
(Beat / Difference Frequency)

We want the Beat freq \approx at 0 so the oscillator has nearly the same freq as precession freq.

1216 Setting X, Y, Z, Z² to zero. Placing mineral oil in RF sample

1223 Obviously something similar to a FID packet (M1) and Difference Signal. See figure 02.11E1




FIGURE 02.11E.1: FID Signal (Channel One) from Oscilloscope with Difference Signal From Q Output (Channnel Two)




FIGURE 02.11E.2: Zoomed in version of the previous figure.




FIGURE 02.11E.3: FID Signal (Channel One) from Oscilloscope with Difference Signals From Q Output (Channnel Two) and Difference Signals From S Output (Channel Three)




FIGURE 02.11E.4: Zoomed in version of the previous figure.

248 Finished pushing images. Taking break. Turning off equipment.

311 Returned, going to do mini commentary on data and then signout/plan experiment for tomorrow

1333 Clearly seeing what seems to be a FID Signal on Channel 1 in Fig 02.11E.1, however it is not a sharp peak as expected from the lab manual [2, pp III-7]. Instead we see a more rounded peak evident in Fig 2.11 E.3. However, Fig 02.11E.4 looks more similar to what is expected. Expected signal in Fig 02.11F.1.





FIGURE 02.11F.1: Expected FID and Difference signals from the Teach Spin Cw Laboratory Manual [2] (p. III-7)

24



~~class discussion break~~

1353 Further investigation and analysis to be conducted

1354 Objectives for next session :

1. Analyze Data From Today
2. Characterize and understand theory

1359

Objectives for lab book :

1. Read and write background theory
2. Read some extractions
3. Plan remaining lab sessions + proposal for Part PS
for overall experiment.

1400

Signing out of lab Reagan

02-14

05 Signifying lab

08

Objectives:

- Discuss Findings w/ P-PS
- Collect data for one more sample*
- Organize (theory should be studied over Reading week).

35

Discussion w/ T-KT

- Data looks good
- Lab manual has interesting extensions
- Lower V_{max} just means pulses will be longer

40

Turning on equipment and trying to obtain same results

45

Parameters changed & from last time. Tuning to Proc (P0211C)

50

Better data obtained!

55

Discussion w/ P-PS

- Overall objectives should be derived
- After Reading week, discuss w/ P-PS about Theory & Ideas.

1058

Adjusted oscilloscope save settings.

1124

Setting putting in light mineral oil.

1127

After cycling scope, light mineral oil data observed. Looks qualitatively different.

1324

Will post images later today. For now signing off at lab ~~for now~~

2025-02-15




FIGURE 02.14A.1: New Data Obtained from scope for heavy mineral oil. We can clearly see a more distinct pulse width forming as is expected from Fig. 02.11F.1 and from the 'spin echo' section of the lab manual [1].



FIGURE 02.14B.1: Peak obtained when measuring light mineral oil. Clearly, as compared to Fig. 02.14A.1, there peak is significantly smaller.

*Time not
recorded*

Significant (just picking images from yesterday)

02-20

Signing in, conducting some light analysis on data collected on 2025-02-20 (pg 26). Mainly just going to write code and verify that things make sense.

20

Writing Analysis code

125

Using normal pd.read_csv('...') for dataframe values.

358

Plots fit very non-ideal, significant issues w/ finding good guess.

701

After significant debugging, there are still issues. Going to discuss w/ Prof PS and/or TA. KT ~~tomorrow~~ on 2025-02-25 (next supervised lab session).

727

Signing out APC

2025-02-23

2157 Signing in

2159 Plan for measuring T_1 and T_2 .22 00 T_2 Plan:

Fit FID signal, this should correspond w/ crude estimates from fitting FSP signal

22 10 From [2] T_1 can accurately be measured using a 2 pulse sequence as

$$\partial_b M_2(t) = (M_0 - M_2(t)) \cdot \frac{1}{T_1}$$

$$\Rightarrow M_2(t) = M_0 \left(1 - \exp \left\{ \frac{-t}{T_1} \right\} \right) \text{ as in the background theory.}$$

$$\pi \xrightarrow{\tau} \pi/2 \xrightarrow{\tau/2} \pi$$

↳ initial amplitude & M_2 pre $\pi/2$ pulse
of FID post
 $\pi/2$ pulse

22 19 So taking a few τ values will give us T_1 .22 20 T_2 can be measured using a reverse sequence $\pi/2 \xrightarrow{\tau} \pi \xrightarrow{\tau} \text{Echo}$ Using $M_m(2\tau) = M_0 \left(\exp \left[-\frac{2\tau}{T_2} \right] \right)$ to determine22 28 echo amplitude as a function of τ yields the T_2 time~~Some graphs~~

22 31 Signing out.

signing into lab

Discussion w/ KT, about future exams

Data Analysis

After debugging, decent initial plot found to tail of pulse.




FIGURE 02.25B.2: Raw Pulse data obtained from the Oscilloscope.
We will castrate this dataset to fit the tail of the exponential.

To fit to this data, castrate $-0.002 \rightarrow 0.000$ and $0.00015 \rightarrow \infty$. This was found through trial and error.

Using this, the values on the following page were found
Page 30

For initial estimation of uncertainties, the min and max of the noise from $-0.002 \leftrightarrow 0.000$ was recorded. This was used as the uncertainty in the y-data. No uncertainty was assumed for x-data.

↑

This will be fine-tuned and analyzed later

Fit to tail of Heavy Mineral Oil Pulse




FIGURE 02.25B.3: Obtained fit using the equation in the graph. Obtained parameters are as follows, $A = 2.0 \pm 70000$, $\tau = -2233 \pm 4$, and $x_0 = 0.0002 \pm 16$ with a $\chi^2_{\text{red}} = 0.2$.

1212

Clearly, we have a decent fit, $\chi^2_{\text{red}} < 1 \Rightarrow$ some overfitting to noise. However, I believe this is much likely an overestimation of uncertainty as the pattern in the residuals is quite insignificant (spacelike far together) and of the pulse width).

1216

It follows then that T_1 for Heavy Mineral Oil is given by

$$\begin{aligned} T_1 &= \frac{1}{\tau} = 0.0004478 \pm 0.0000008 \\ &= 447.8 \pm 0.8 \text{ ms} \end{aligned}$$

1221

Same procedure will be followed for Light Mineral Oil.

1222

Sighing out at Vb again

2-28

Sizing into lab

Discussion w/ LT-Carry about old NMR machine
- Issue probably due to shielding on copper tip of probe
becoming desoldered

0

(Previous) Discussion w/ Arvish (prior student doing NMR).
- Observed problem w/ old NMR machine is that no Spin Echo
was observed.

25

35 Discussion w/ TA-KT:

- Should try and observe spin echo
- No expected result if we don't sample w/ laser
- Fixing old machine would be interesting.

46

Helped another setup briefly.

053

Lab Objectives:

- Measure T_2 ,
- Try and observe Spin Echo
- Setup Spin Diffusion in ~~H₂O~~ Distilled Water.

055

T_2 Measurement Procedure (P02.28A)

↑ π pulse

1. Tune Spectrometer to Resonance and obtain best π pulse possible.

2. Adjust Fine Tuning Capacitor, if FID tail smaller => keep adjusting in that direction.

3. Adjust till no tail (or nearly no tail) observed in FID signal.

1105

Scope settings loaded from 2025-02-14, no signal observed.
Checking apparatus parameters from P02.11C

1116

New Signal observed, having different from before.




FIGURE 02.28A.1: Different signal than previous measurements observed. Need to retune and adjust parameters.

1130 Break

1146 Returned.

Adjusting to obtain π -pulse parameters for Heavy Mineral oil as per {PO2.28A}.

1148 Discussion w/ Prof. Tahir Shaaran

{ St=95% of gradient (plot $\ln(I/I_0)$ against G^*)

and use Stejskal-Danner Equation to obtain spin diffusion coefficient

→ For measuring spin diffusion in distilled water.

1152 Returning to π -pulse parameters [11.48.2025-02-28]
1157 Reading TechSpin manual.

1210 B-pulse not working, procedure outlined in manual not working.

1211 Adjusting ~~manually~~ w/o manual

1236 Spin Echo ● Observed?




FIGURE 02.28B.1: Possible Spin Echo observed. First peak is from a $\frac{\pi}{2}$ -pulse, second peak is a π -pulse. The third peak may be the observation of a spin echo.

1240 Parameters (PP):

A_len : 3.90

Sync: A

B_len : 7.80

A_t: on

T : 0.0042

B_t: on

Num_B: 1

M_t: off

Period: 100.0 ms

Rest all same as before.

1241 Having a headache. Going to take a break and do reading + Analysis over weekend.

304 signing off to bed Afri-

2025-03-03

1514

Signing into lab

Measure T_2^*

1520 discussion w/ Shrawan K. Shridhar (SKS)

Testing = T then T_2 pulse.

st

1528 Scope 45 is hist^{Image.} for today

→ SCROFULOUS PULSES
↳ third order invariant

1535 SKS suggested FFT give you beat frequency allowing us to find resonance.

1547 Fixing FFT code.

1604 Nell to All in $f: \approx 20.86200 \text{ MHz}$

1605 Survival Br.

03-07)

Siging in

Explanation of work conducted today @ 1520:

- Discussion w/ SKS about data
- SKS noticed pulse seemed to be off Resonance and suggested a method to find Resonance

B:23

Take Q or I data and get Fourier Transform. Out put should be similar to Fig 03.03A.1




FIGURE 03.03A.1: Data obtained from taking the fourier transform of the I signal data. Notice the existance of two distinct peaks. One being the mean and the other being the resonance frequency.

837

Now take ~~Maxima~~ mean and remove it from the signal. The next highest peak should be near the resonance frequency. This can be seen in Fig 03.03A.2. ~~The~~




FIGURE 03.03A.2: Mean removed fourier transform.

18:42

Finding this resonance frequency has the coil and applying it to a
an edit of f_R to the ~~Resonance~~ machine has the following shown on
the I pulse (Fig 03.03B.1)




FIGURE 03.03B.1: Resultant I-signal observed after tuning using fourier transofmrr method suggested by SKS.

18:52

Further edits/pine adjustments yielded Fig 03.03B.2. From here,
 $f_R = 20.8620$ MHz was found as the Resonance frequency at
Heavy Mineral oil.

19:00

Lab Objectives for 2025-03-04

1. Rebutte HMO and LMO data.
2. Record T_1 and $\textcircled{T}_2 \text{ times}$ using π then $\pi/2$ pulse
Spin Relaxation
3. Look into spin relaxation for mixing operators?

~~Typing out & ref~~

Note, these notes were filled in lab. Images were posted
the following week day





FIGURE 03.03B.2: Corresponding signal.

1905 Signing out APR



S-03-04

30

Sign in to lab

131

Tuning on equipment, reviewing notes from 2025-03-01

948

~~Finding~~ Tuning Equipment

952

Tuned to $f_R = 20.88620 \text{ MHz}$, $\Delta f_{\text{min}} = \cancel{\Delta f} = \cancel{\Delta f}$
~~Had to perform fine adjustments~~

$$\ell_{\pi_1} = 3.90 \mu\text{s}, \ell_{\pi_2} = 7.80 \mu\text{s}$$

0959

Lab Objectives

1. Verify f_R
2. Release data for HMO and LMO
3. Take data for distilled H2O

1013


$$\text{Verified } \ell_{\pi_1} = 3.80$$

Difference

$$\Rightarrow \ell_{\pi} \approx 7.60$$

→ we may verify this now

1014

Output w/ $\ell_{\pi} = 6.94 \mu\text{s}$ in Fig. 03.04A.1FIGURE 03.04A.1: Result of a π pulse with pulse length $\ell_{\pi} = 6.94 \mu\text{s}$.

016 Noticing there is still an unexpected FID signal. Consulting lab manual [1] to verify procedure.

027 New parameters (which I hope are correct):

$$t_{\pi_0} = 3.80 \mu s \pm 10 \mu s$$

$$t_\pi = 7.40 \mu s \pm 10 \mu s$$

$$f_R = 20.86114 \text{ MHz}$$

Relaxant output in Fig 03.04B.1




FIGURE 03.04B.1: Relaxant output with new pulse length parameters.

030 Now, we may attempt to observe the spin distillation in water. Relaxant theory will be added later. For now, important immediate information follows:




Fig 03.04C.1




Fig 03.04C.2




Fig 03.04C.3

Fig 03.04C: Bloch spheres of ground state (C1), $\pi/2$ pulse (C2) and π pulse (C3)

1033

In order to measure the spin distribution in sp₃ distilled H₂O, consider the following procedure:

1. Flip spins from Ground state,
2. Interrogate Spins using $\pi/2$ pulse.
3. Measure Relaxation time.

1050

Sample filled ~~out~~ w/ distilled H₂O. Putting on
Machine.

1055

Sample stuck in Machine, asking LT for help

1058

Was able to remove sample. Setting machine up again

1101

Relaxing output signal very small. Need to feed in signal into op-amp

1102

Coffee Break (25 min max) — *Did I end up taking*



FIGURE 03.04D.1: Output signal from distilled water sample. Notice that the signal strength is quite low and very noisy.

just stabilized amplification device is a lock-in amplifier, not a non-inverting op-amp.

1112 trying to quickly make non-inverting op-amp using LM358




Fig 03.04 E1

$$G_{\text{gain}} = 1 + \frac{R_1}{R_2}$$

↑ closed loop voltage gain,

We want ≈ 10 gain

$$\Rightarrow R_1 = 1 \text{ k}\Omega,$$

$$R_2 = 100 \Omega .$$





Fig 03.04 E2

LM358 Pinout




Fig 03.04 E3

- | | |
|---------------------------|--|
| $1 = X$ | |
| $2 = V - \text{in}$ | |
| $3 = V + \text{in}$ | |
| $4 = X$ | |
| $5 = X$ | |
| $6 = \text{Power (+) in}$ | |
| $7 = V_{\text{out}}$ | |
| $8 = V_{\text{out}}$ | |

1147 Designed and finished soldering op-amp. Testing break and shorting link to Pin 4. A.T.H.

1219

OpAmp Remade




FIGURE 03.04F.1: Breadboard of non-inverting op-amp.

1220

Still not fully working as intended, taking $\text{d}m$ into
and going to delay circuit at home.

1220

Setting pulse delay $\tau = 50\mu\text{s}$, Running π followed by $T/2$ pulse
and recording d_m .

1222

Signal extremely noisy




FIGURE 03.04G.1: Noisy small signal observed.

1222 Adjusting delay, switch Escape to time mode averaging.

1224 Results not good, switching to High-Res Mode.

1225 Resulting signal observed




FIGURE 03.04G.2: Observation for Distilled Water.

1226 Adjusting $\tau = 5\mu s$, so that $\pi/2$ pulse $T = 7.85 \mu s \approx 8\mu s$
so that $\pi/2$ pulse directly follows π pulse. Resulting signal




FIGURE 03.04G.3: Observation for Distilled Water.

1230

Turning Blanking off has the following output




FIGURE 03.04G.4: Observation for Distilled Water.

1231

Not conducting Data Analysis Now/Today. Will now rebake some Heavy & Light Mineral Oil Measurements after a break

1232

Taking Break.

1255

Returned from break

1255

Taking Readings for HMO and LMO

Heavy Mineral Oil
Light Mineral Oil

1256

Removed Distilled H₂O, put in LMO

1257

Gain reduced to 0.75 as before.




FIGURE 03.04H.1: Observation for Light Mineral Oil.

1258 ~~PS~~ Changing h π pulse




FIGURE 03.04H.2: Observation for Light Mineral Oil.

1259 π pulse followed by $\frac{2\pi}{2}$ pulse w/ $\tau = 0.0007$ s




FIGURE 03.04H.3: Observation for Light Mineral Oil.

1300 Presumably, Spin Echo observed.

1301 Changing sample to Heavy Mineral Oil and repeating.

1302 Heavy Mineral Oil inside machine.




FIGURE 03.04I.1: Observation for Heavy Mineral Oil.

B02 π -pulse




FIGURE 03.04I.2: Observation for Heavy Mineral Oil.

B03 π followed by $\tau\pi_2$ pulse




FIGURE 03.04I.3: Observation for Heavy Mineral Oil.

303 Data collected. Need to process images and paste them into book.

303 Objectives for today/tomorrow:

1. Finish writing Background Theory
2. Write up lit review on applications
3. Data Analysis

1304 May need to email supervisor for an extension, we'll see.

1108 Checking lab manual incase any details were missed

1309 There is some interesting discussion on using the PS2 controller to "keep the spins alive".

1310 Let's see if we can keep the spins alive and do some sort of unitary operation like in a G Algorithm.

311 Maximum measurement at $t = 0.0004$ s




FIGURE 03.04J.1: Playing around with pulses to see if we can implement a unitary gate. This is more just exploration with no expectation of results.

1318

Running only $\frac{\pi}{2}$ pulse and seeing if we can extract FID.

1324

Interesting effects observed, but not much worth reporting.

1325


Main interesting change was that, apparently, the magnetic field gradient in the $+z$ direction had some small value. After changing this f_z changed to $f_z = 20.85558 \text{ MHz}$. However, due to time constraints, we will still analyze old data.

1326

We may do this because that data was still taken consistently and any comparative analysis will still be self-consistent.

1327

Back to exploring. Consider the following algorithm/circuit



329

where \hat{H} is the Hadamard operator w/ unitary matrix

$$\hat{H} = \frac{1}{\sqrt{2}} \begin{bmatrix} 1 & 1 \\ 1 & -1 \end{bmatrix}$$

and \oplus is a "CNOT" or controlled NOT operator. This may also be ~~written~~ called a controlled pauli X operator

$$\hat{X} = \begin{bmatrix} 0 & 1 \\ 1 & 0 \end{bmatrix}$$

1331 We may examine this circuit w/ initial ket at $|1\rangle_0$


$$|\Psi_0\rangle = \cancel{|100\rangle}$$

For \hat{H} acting on Q_0

$$\begin{aligned}\hat{H}|0\rangle &\mapsto \frac{1}{\sqrt{2}}(|0\rangle + |1\rangle) \\ \hat{H}|1\rangle &\mapsto \frac{1}{\sqrt{2}}(|0\rangle - |1\rangle)\end{aligned}$$

$$\text{so } |\Psi_1\rangle = \frac{|100\rangle + |01\rangle}{\sqrt{2}}$$

1347 Now the CNOT gate works as follows




$$\begin{aligned}\text{Then } |\Psi_2\rangle &= \text{CNOT}(|\Psi_1\rangle) \\ &= \frac{|100\rangle + |111\rangle}{\sqrt{2}}\end{aligned}$$

Notice this is ρ_B ie, a Bell state or an EPR pair.

349

We obviously cannot do such manipulations in ~~the~~ lab, however, we can model a CNOT gate as follows.

If we model our groundstate as pre π-pulse, then π-pulse = \hat{X} gate



1351

I would like to see if we can take a measurement and if the state was $|0\rangle$ or $|1\rangle$.

1352

An ideal π -pulse has no FID, but this is not ideal. If we interrogate the spin w/ a $\frac{\pi}{2}$ pulse, we may see see a difference b/w a $|0\rangle$ & $|1\rangle$.

1354

Data obtained looks normal (spin echo).
~~takes long~~ $\pi - \tau - \pi_{1/2} - \boxed{\pi}$

We may also try $\pi - \tau - \pi - \pi_{1/2} - \boxed{\pi}$

1358

Nothing observed (need a better way to imprint).

1359

Maybe we can try a H gate

$$\hat{H} = \frac{\hat{X} + \hat{Z}}{\sqrt{2}} \quad \hat{X}\hat{Z} = \begin{bmatrix} 0 & 1 \\ 1 & 0 \end{bmatrix} \begin{bmatrix} 0 & 1 \\ 1 & 0 \end{bmatrix} = \boxed{0}$$

$$\hat{E} X = \hat{H} \hat{Z} \hat{H}$$

Not today ~~as~~, ~~long~~ will think more, for now need to get back on track.

1400

Done for today, turning off equipment. taking a 2nd group circuit to debug.

1402

Signing out of lab ~~for now~~


- 1920 Signing in.
- 1923 ~~↳ Objectives:~~
1. Discuss Applications of NMR from a few papers.
 - ↳ NMR (MRI) in Medicine
 - ↳ Quantum Computing.
- 1927 Given my background, I will first start w/ a review for Quantum (a) Computing.
- 1927 From [6] it was reported that "143 is the largest number yet factored by a Quantum Algorithm." What does this actually mean and what does it have to do w/ NMR?
- 1937 In general, the problem of factoring is not so glorious, it is indeed one of the first things we're taught in any schooling. Factoring 36 for instance is not so impressive: $(1, 36)$, $(2, 18)$, $(3, 12)$, $(4, 9)$ and $(6, 6)$. However, most of modern encryption is built on the idea that it is really really hard to factor prime numbers. For composite numbers (like 36) we may use techniques such as the Sieve or Euclidean or Probabilistic Factoring Algorithms to achieve factorizations. It is also not so difficult to check if a number is a dead prime. Wilson's Algorithm, for instance, is a deterministic test that can verify a number is prime in $O(n^2)$ time complexity.
- Now
- However, as the size of the numbers scale, it becomes increasingly complex to factorize these numbers. Encryption methods (for instance RSA or SHA256) take two large primes and ~~factorize~~ multiply them to obtain a composite number with only two factors (other than one and itself), which are two primes.
- 1920 "When an integer reaches more than 100,000 or so digits, it becomes too complex to ~~factorize~~ ~~classically~~ solve using classical computing methods." [6]. Quantum Computers on the other hand have some "slick tricks" which they can use to overcome this hurdle.
- 1923 Shor's Algorithm, probably the most famous quantum algorithm (apart from maybe Grover's Algorithm) utilizes the speed advantage exclusive to quantum computing like entanglement to leverage the "Quantum Advantage!"

→
2024 Shor's Algorithm effectively works as follows:

Let $N = ab$, select c such that c^{-1} is coprime to $M = (a-1)(b-1)$.
Now, compute d such that $c \cdot d \bmod M = 1$.

Choose ~~odd~~ d which is coprime to N ($\gcd(d, N) = 1$). Then $f_{N,d}$ which is periodic with period p . Then, we may find x if
Either $\gcd(d^{p_1} + 1, N)$ or $\gcd(d^{p_2} - 1, N)$ is a factor of N .

2024 And the corresponding circuit diagram [7]:



2024 I will not be going in depth into how the Ansatz or Oracle work in this case.

2025 It was found, as highlighted by the article, that in 2001, E. Farhi, Fetter proposed an algorithm known as AQC or Adiabatic Quantum Computing.


2027 Unlike Shor's Algorithm, AQC relies on quantum adiabatic processes which all possible solutions to the factoring problem are encoded as Eigenstates to our problem Hamiltonian.

2111 In the physical implementation, ~~gradually~~ NMR is used to gradually evolve the Hamiltonian to solve the factoring problem.


2137 I realize now that this article does not go nearly in-depth enough to explain how NMR works in this context, so I will attempt to supplement.

2140 One implementation of gate operations utilizes the energy ~~of~~ spin states to encode one qubit implementation.

2141 Consider the following implementation of a ~~control~~ CNOT gate.



$$\begin{aligned}
 |00\rangle &\rightarrow |00\rangle \longrightarrow |00\rangle \longrightarrow |00\rangle \\
 |01\rangle &\rightarrow |01\rangle \longrightarrow |01\rangle \longrightarrow |01\rangle \\
 |10\rangle &\longrightarrow |00; m\rangle \longrightarrow -|00; m\rangle \longrightarrow |10\rangle \\
 |11\rangle &\longrightarrow |01; m\rangle \longrightarrow |01; m\rangle \longrightarrow -|11\rangle
 \end{aligned}$$



Note the above is ~~for~~ for a Trapped Ion system. Notice the use of π and 2π pulses to manipulate b/w energy states and moving states (denoted by the m). In this way we can use NMR - esqun techniques to implement 2-qubit gates and other unitary operations.

2003 signing act. ~~ATC~~

2025-03-05

14/17

Signing in


14/19

Continuing Analysis of collected data on LMO, HMO, and d-H₂O from yesterday (2025-03-04).

14/25

I have realized an error in my data collection, I only took one T for each measurement, making it impossible to properly fit T_1 and T_2 . I can still get an estimate using my FID code from earlier perhaps.

Fit to tail of lmo Pulse



Residuals






FIGURE 03.05
03.05
03.05
FIGURE 03.05A.1: $A = 125.02730736861422 \pm \text{inf}$, $\tau = -159.21440895422927 \pm \text{inf}$, $x_0 = -0.711144259571313 \pm \text{inf}$ $\chi^2_{\text{red}} = 2.0$

14/25

Note that due to personal time constraints most of these plots are not optimized. Especially Fig 03.05C.1 seems to be fit to the wrong data.



03.05
FIGURE 03.05B.1: $A = 0.0 \pm 9000.0$, $\tau = -5456.4 \pm 90.0$, $x_0 = 0.0 \pm 200.0$, $\chi^2_{\text{red}} = 0.003$



03.05
FIGURE 03.05C.1: $A = 0.650974392642326 \pm 1436422.7844434683$,
 $\tau = 2375.446650108936 \pm 224.61004122688507$, $x_0 = 0.0003897026532702442 \pm 937.2177888340759$, $\chi^2_{\text{red}} = 2.0$

1510

Finally I will briefly discuss what I should have done.
 Namely I may, set 5C periods and varied them

1511

This would have given me sufficient data to properly calculate T_1 and T_2 as currently my datasets, and ~~had time~~ are too restricted to get meaningful results.

1513

I will comment that the uncertainties in Fig 03.05A.1, Fig 03.05B.2, and Fig 03.05C.1 are quite unrealistic. Here my methodology was simply based off of the noise and taking the min and max of those noisy ~~percent~~ areas.

1517

Likewise, taking ~~and~~ another FFT and taking the high order parameters may have been a better method for such analysis of uncertainties.

1521

One primary error was ~~not~~ realizing I had ~~an~~ ~~an~~ incorrect resonance condition early on. This lead to rushed final data collection, ~~and~~ ~~it~~ should have been more time, I would have realized my mistake sooner and taken more data.

```

#!/usr/bin/env python3
# -*- coding: utf-8 -*-
"""
Created on Tue Feb 25 10:25:15 2025

@author: Aditya K. Rao
@github: @adirao-projects
"""

import numpy as np
import pandas as pd
import matplotlib.pyplot as plt
import toolkit as tk
import toolkitnew as tkn

def load_data(csvname):
    # ch1 is ENV out
    # ch2 is Q out
    # ch3 is I out

    df = pd.read_csv(csvname, names=['Time', 'ENV', 'Q', 'I'], skiprows=2)
    df_units = {'Time': 's', 'ENV': 'V', 'Q': 'V', 'I': 'V'}
    return df, df_units

def pulse(x, y0, x0, w0, h, tau):
    if x <= x0:
        return 0

    elif x <= x0 + w0:
        return h

    else:
        return h * np.exp((x - x0) * tau) + y0

def fiteqn(t, M0, T1):
    return M0 * (1 - np.exp(-t / T1))

def get_uncert(df):
    df_ = df.where(df['Time'] < 0)
    df_['uENV'] = df_[['ENV']].max() - df_[['ENV']].min()
    df_['uQ'] = df_[['Q']].var()
    df_['uI'] = df_[['I']].var()

    return df

def fidcalc(df, units, name):
    df = df[df['Time'] < 0.001]
    df = df[df['Time'] >= 2e-4]

    meta = {'title': f'Fit to tail of {name} Pulse',
            'xlabel': 'Time (s)',
            'ylabel': 'ENV (V)',
            'fit-label': r"$A \exp(-\tau \ln(x - x_0))$",
            'data-label': "Pulse Data",
            'loc': 'upper right',
            'save-name': f'{name}fid'}
```

```

tkn.quick_analyze(df['Time'], df['ENV'], yerr=df['uENV'], fit_type='exp',
    res=True, chi=True, guess=(1.75, -100, 0.00025),
    meta = meta, show=True, save=True, dataout=True,
    params=('A', r'\tau', 'x_0'), rounding=False)

def t1calc(df, units, name):
    df = df['Time']<0.001
    df = df[df['Time']>=2e-4]
    vpulse = np.vectorize(pulse)

    meta = {'title':f'Fit to tail of {name} Pulse',
            'xlabel':'Time (s)',
            'ylabel':'ENV Voltage (V)',
            'fit-label': r"$M_0(1 - \exp(-\frac{t}{T_1}))$",
            'data-label': "Pulse Data",
            'loc':'upper right',
            'save-name':f'{name}.png'}

    tkn.quick_analyze(df['Time'], df['ENV'], xerr=df['uENV'], fit_type='custom',
        model_function_custom = fiteqn, res=True, chi=True,
        meta = meta, show=True, save=True, dataout=True,
        params=('A', r'\tau', 'x_0'))


def t2calc(df, units, name):
    df = df[df['Time']<0.001]
    df = df[df['Time']>=2e-4]
    vpulse = np.vectorize(pulse)

    meta = {'title':f'Fit to tail of {name} Pulse',
            'xlabel':'Time (s)',
            'ylabel':'ENV (V)',
            'fit-label': r"$A\exp(-\tau(x-x_0))$",
            'data-label': "Pulse Data",
            'loc': 'upper right',
            'save-name':f'{name}'}

    tkn.quick_analyze(df['Time'], df['ENV'], xerr=df['uENV'], fit_type='exp',
        res=True, chi=True, guess=(1.75, -100, 0.00025),
        meta = meta, show=True)

if __name__ == '__main__':
    hmo = ['./Data/hmo1.csv', './Data/hmo2.csv', './Data/hmo3.csv']
    lmo = ['./Data/lmo1.csv', './Data/lmo2.csv', './Data/lmo3.csv']
    h2o = ['./Data/h2o1.csv', './Data/h2o2.csv', './Data/h2o3.csv']

    for sample_name, sample_dat in zip(hmo, lmo, h2o):
        for data_type, data_path in zip('pi2', 'pi', 'pipi2'), sample_dat):
            df, units = load_data(data_path)
            df = get_uncert(df)
            df = df.dropna()

            if data_type == 'pipi2':
                pass

```

```
#tIcalc(df, units, sample_nam)  
elif data_type == 'pi2':  
    fidcalc(df, units, sample_nam)
```

```
# -*- coding: utf-8 -*-
"""
Created on Tue Jan 23 12:34:34 2024
Updated on Mon Oct 14 21:22:09 2024
Updated on Mon Feb 10 09:51:03 2025
Updated on Sat Mar 01 21:03:54 2025
Updated on Mon Mar 03 14:45:13 2025

```

Lab Toolkit

@author: Aditya K. Rao
 @github: @adirao-projects
 one

```
import numpy as np
from scipy.optimize import curve_fit
import matplotlib.pyplot as plt
import matplotlib.gridspec as gridspec
import os
import math
import textwrap

#from uncertainties import ufloat

font = {'family': 'DejaVu Sans',
        'size' : 30}

plt.rc('font', **font)

def fit_data(xdata, ydata, fit_type, override=False,
             override_params=(None,), uncertainty=None,
             res=False, chi=False, uncertainty_x=None,
             model_function_custom=None, guess=None):

    def chi_sq_red(measured_data:list[float], expected_data:list[float],
                   uncertainty:list[float], v: int):
        if type(uncertainty)==float:
            uncertainty = [uncertainty]*len(measured_data)
        chi_sq = 0

        # Converting summation in equation into a for loop
        # for i in range(0, len(measured_data)):
        #     chi_sq += (pow((measured_data[i] \
        #                     - expected_data[i]),2)/(uncertainty[i]**2))

        for m, e, u in zip(measured_data, expected_data, uncertainty):
            chi_sq += ((m-e)/u)**2

        chi_sq = (1/v)*chi_sq

        return chi_sq

    def residual_calculation(y_data: list, exp_y_data) -> list[float]:
        residuals = []
        for v, u in zip(y_data, exp_y_data):
            residuals.append(u-v)

        return residuals
```

```

def model_function_linear_int(x, m, c):
    return m*x+c

def model_function_exp(x, a, b, c):
    return a*np.exp(b*(x-c))

def model_function_log(x, a, b):
    return b*np.log(x+a)

def model_function_linear_int_mod(x, m, c):
    return m*(x+c)

def model_function_linear(x, m):
    return m*x

def model_function_xlnx(x, a, b, c):
    return b*x*(np.log(x)) + c

def model_function_ln(x, a, b, c):
    return b*(np.log(x)) + c

def model_function_sqrt(x, a):
    return a*np.sqrt(x)

model_functions = [
    'linear' : model_function_linear,
    'linear-int' : model_function_linear_int,
    'xlnx' : model_function_xlnx,
    'log' : model_function_log,
    'exp' : model_function_exp,
    'custom' : model_function_custom
]

try:
    model_func = model_functions[fit_type]
except:
    raise ValueError(f'Unsupported fit-type: {fit_type}')

if not override:
    new_xdata = np.linspace(min(xdata), max(xdata), num=100)

    if type(uncertainty) == int:
        abs_sig = True
    else:
        abs_sig = False

    if guess is not None:
        popt, pcov = curve_fit(model_func, xdata, ydata, sigma=uncertainty,
                               maxfev=20000, absolute_sigma=abs_sig, p0=guess)
    else:
        popt, pcov = curve_fit(model_func, xdata, ydata, sigma=uncertainty,
                               maxfev=20000, absolute_sigma=abs_sig)
    param_num = len(popt)

```

```

meta = {'title' : 'INSERT-TITLE',
        'xlabel' : 'INSERT-XLABEL',
        'ylabel' : 'INSERT-YLABEL',
        'chisq' : 0,
        'fit-label': "Best Fit",
        'data-label': "Data",
        'save-name' : 'IMAGE',
        'loc' : 'lower right'}

main_fig.set_title(meta['title'], fontsize = 46)
if len(uncertainty_x)==0:
    main_fig.errorbar(xdata, ydata, yerr=uncertainty, #xerr=uncertainty_x,
                       markersize='4', fmt='o', color='red',
                       label=meta['data-label'], ecolor='black')
else:
    main_fig.errorbar(xdata, ydata, yerr=uncertainty, xerr=uncertainty_x,
                       markersize='4', fmt='o', color='red',
                       label=meta['data-label'], ecolor='black')

main_fig.plot(plot_x, plot_y, linestyle='dashed',
              label=meta['fit-label'])

main_fig.set_xlabel(meta['xlabel'])
main_fig.set_ylabel(meta['ylabel'])
main_fig.legend(loc=meta['loc'])

if res:
    res_fig.errorbar(xdata, residuals, markersize='3', color='red', fmt='o',
                      yerr=uncertainty, ecolor='black', alpha=0.7)
    res_fig.axhline(y=0, linestyle='dashed', color='blue')
    res_fig.set_title('Residuals')

if save:
    plt.savefig(f"figures/{meta['save-name']}")

if show:
    plt.show()

def block_print(data: list[str], title: str, delimiter='=') -> None:
    """
    Prints a formated block of text with a title and delimiter
    Parameters
    -----
    data : list[str]
        Text to be printed (should be input as one block of text).
    title : str
        Title of the data being output.
    delimiter : str, optional
        Delimiter to be used. The default is '='.
    Returns
    -----
    None.
    Examples
    -----
    >>> r.log = 100112.24998716181

```

```

>>> r_dec = 0.007422298127465114
>>> data = [f'r^2 value (log): {r_log}', 
           f'r^2 value (real): {r_dec}']
>>> block_print(data, 'Regression Coefficient', '=')
===== Regression Coefficient =====
r^2 value (log): 100114.24998718781
r^2 value (real): 0.007422298127465114
=====
"""
term_size = os.get_terminal_size().columns

breaks = 1
str_len = len(title)+2
while str_len >= term_size:
    breaks += 1
    str_len = math.ceil(str_len/2)

str_chunk_len = math.ceil(len(title)/breaks)
str_chunks = textwrap.wrap(title, str_chunk_len)
output = ''
for chunk in str_chunks:
    border = delimiter*(math.floor((term_size - str_chunk_len)/2)-1)
    output = f'{border} {chunk} {border}\n'

output=output[:-1]

output+= '\n' + '\n'.join(data) + '\n'
output+=delimiter*term_size

print(output)

def quick_analyze(xdata, ydata, fit_type, model_function_custom=None,
                  yerr=None, xerr=[], res=False, chi=False, guess=None, \
                  save=False, meta=None, dataout=False, params=None,
                  rounding=True, show=False):

    data = fit_data(xdata, ydata, fit_type=fit_type, uncertainty=yerr,
                    chi=chi, res=res, guess=guess,
                    model_function_custom=model_function_custom)

    plot_data(xdata, ydata, data['plotx'], data['ploty'],
              residuals=data['residuals'], uncertainty=yerr,
              uncertainty_x=xerr, res=res, save=save, meta=meta,
              show=show)

    if dataout:
        with open(f'figures/{meta["save-name"]}.dat', 'w') as f:
            for n, t, d in zip(params, data['nopt'], data['pstd']):
                if rounding:
                    d = round(d, -int(np.floor(np.log10(np.abs(d))))) 
                    dp = len(str(d).split('.')[1])
                    t = round(t, dp)

                f.write(f'{n} : ${n} = {t} \pm {d}$\n')
            chisq = round(data['chisq'], \

```

```

    -int(np.floor(np.log10(np.abs(data['chisq'])))))
f.write(r'chisq : $\chi^2_{\text{red}}$ = ' + str(chisq) + '\n')

def param_print(names, popt, pstd, latex=True):
    """
    Prints all parameters for easy input into a LaTeX document.

    Parameters
    -----
    names : list[str]
        Symbols/Names of each parameter
    popt : list[float]
        Value of each parameter.
    pstd : list[float]
        Standard deviation/uncertainty associated with each parameter.
    latex : bool, optional
        Print in LaTeX format. The default is True.

    Returns
    -----
    None.
    """

    data = []
    for n, t, d in zip(names, popt, pstd):
        if latex:
            data.append(f'{n} = {t} \pm {d}')
        else:
            data.append(f'{n} = {t} +/- {d}')
    block_print(title='Parameters', data=data)

def numerical_methods(method_type, args=None, custom_method=None):
    def gaussxw(N):
        """
        Initial approximation to roots of the Legendre polynomial
        x = np.linspace(1,4*N-3,N)/(4*N+2)
        x = np.cos(np.pi*a+1/(4*N*N*np.tan(a)))

        Find roots using Newton's method
        epsilon = 3e-16
        delta = 1.0
        while delta>epsilon:
            p0 = np.ones(N,float)
            p1 = np.copy(x)
            for k in range(1,N):
                p0, p1 = p1, ((1-x*k)*x*p1-k*x*p0)/(k+1)
            dp = (N+1)*(p0-x*p1)/(1-x*x)
            dx = p1/dp
            x -= dx
            delta = max(abs(dx))
        # Calculate the weights
        w = 2*(N+1)*(N+1)/(N*x*x*(dp*dp))
        """

```

```
    return x, w

def gaussxwab(N,a,b):
    x,w = gaussxw(N)
    return 0.5*(b-a)*x+0.5*(b+a),0.5*(b-a)*w

methods = {
    'gaussxw' : gaussxw,
    'gaussxwab' : gaussxwab,
    'custom' : custom_method
}

try:
    method = methods[method_type]

except:
    raise ValueError(f'Unsupported method-type: {method_type}')

return method(*args)

def interpolation_methods(method_type, args=None, custom_method=None):

methods = {
    'custom' : custom_method
}

try:
    method = methods[method_type]

except:
    raise ValueError(f'Unsupported method-type: {method_type}')

return method(*args)
```

```

#!/usr/bin/env python3
# -*- coding: utf-8 -*-
"""
Created on Tue Feb 25 10:25:15 2025
@author: Aditya K. Rao
@github: @adirao-projects
"""

import numpy as np
import pandas as pd
import matplotlib.pyplot as plt
import toolkit as tk
import toolkitnew as tkm

def load_data(csvname):
    # ch1 is ENV out
    # ch2 is Q out
    # ch3 is I out

    df = pd.read_csv(csvname, names=['Time', 'ENV', 'Q', 'I'], skiprows=2)
    df_units = {'Time': 's', 'ENV': 'V', 'Q': 'V', 'I': 'V'}
    return df, df_units

def pulse(x, y0, x0, w0, h, tau):
    if x <= x0:
        return 0

    elif x <= x0 + w0:
        return h

    else:
        return h * np.exp((x - x0) * tau) + y0

def get uncert(df):
    df_ = df.where(df['Time'] < 0)
    df['uENV'] = df_['ENV'].max() - df_['ENV'].min()
    df['uQ'] = df_['Q'].var()
    df['uI'] = df_['I'].var()

    return df

def get_t1(df, units, name):
    df = df[df['Time'] < 0.001]
    df = df[df['Time'] >= 2e-4]
    vpulse = np.vectorize(pulse)

    # data = tk.curve_fit_data(df['Time'], df['ENV'], uncertainty=df['uENV'],
    # fit_type='custom', chi=False, res=True,
    # model_function_custom=vpulse,
    # guess=(0.95, 0.00505, 0.0062, 1.75, -0.0001))
    data = tk.curve_fit_data(df['Time'], df['ENV'], uncertainty=df['uENV'],
    fit_type='exp', chi=True, res=True,
    guess=(1.75, -100, 0.00025))

    print(data['pept'])

```

```

meta = {'title':f'Fit to tail of {name} Pulse',
        'xlabel':'Time (s)',
        'ylabel':'ENV (V)',
        'fit-label': r"$A\exp(-\tau(x-x_0))$",
        'data-label': "Pulse Data",
        'loc': 'upper right',
        'save-name':f'{name}.png'}

# tk.block_print(data=[f"Chi:{data['chisq']}\nPM", f"A:{data['popt'][0]}",
#                   f"B:{data['popt'][1]}", f"x_0:{data['popt'][2]}"],
#                   title='Fit Parameters')
tk.parameter_print(['A', r'\tau', 'x_0'], data['popt'], data['pstd'])

print(r'\chi^2_{\nu} = '+str(data['chisq']))
print(f'T_1 = {np.abs(1/data['popt'][1])}')
tk.quick_plot_residuals(df['Time'], df['ENV'],
                        data['plotx'], data['ploty'],
                        uncertainty=df['uENV'],
                        residuals=data['residuals'],
                        meta=meta)

def get_tinew(df, units, name):
    df = df[df['Time']<0.001]
    df = df[df['Time']>=2e-4]
    vpulse = np.vectorize(pulse)

    meta = {'title':f'Fit to tail of {name} Pulse',
            'xlabel': 'Time (s)',
            'ylabel': 'ENV (V)',
            'fit-label': r"$A\exp(-\tau(x-x_0))$",
            'data-label': "Pulse Data",
            'loc': 'upper right',
            'save-name':f'{name}.png'}

    tkn.quick_analyze(df['Time'], df['ENV'], xerr=df['uENV'], fit_type='exp',
                      res=True, chi=True, guess=(1.75, -100, 0.00025),
                      meta = meta)

def get_fft(df):
    print('here')

    print(df['I'].head())
    i_dat = df['I'].to_numpy()

    #df.plot(x='Time', y='Q')
    i_fft = np.fft.rfft(i_dat)
    res = df['Time'].iloc[1] - df['Time'].iloc[0]
    #print(q_fft)
    #print(res)

    print(i_dat.shape[0])
    freq = np.fft.rfftfreq(i_dat.shape[0], res)

```

```

    freq = np.arange(0, i_dat.shape[0], res)

    plt.plot(freq, np.abs(i_fft), color='black')

    i_dat_new = i_dat - np.mean(i_dat)

    i_fft_new = np.fft.rfft(i_dat_new)

    i_fft_new_abs = np.abs(i_fft_new)

    plt.plot(freq, i_fft_new_abs, color='black', linestyle='dashed')

    zimax = np.max(i_fft_new_abs)

    fmax = np.argmax(i_fft_new_abs)

    print(freq[fmax])

    print(20.9e6 + freq[fmax])
    print(20.9e6 - freq[fmax])

    plt.show()

if __name__ == '__main__':
    df_hmol, units_hmol = load_data('../Data/heavyminoil.csv')
    #df_hmol.plot(x='Time', y='ENV', legend=False, grid='on',
    #              title='Raw Pulse Data for Heavy Mineral Oil',
    #              xlabel='Time (s)', ylabel='ENV (V)', figsize=(14, 8))
    df_hmol = get_uncert(df_hmol)

    #get_fft(df_hmol)

    #print(df_hmol.head())

    #t1_hmolnew = get_t1new(df_hmol, units_hmol, 'Heavy Mineral Oil')
    t1_hmol = get_t1(df_hmol, units_hmol, 'Heavy Mineral Oil')

    #df_lmol, units_lmol = load_data('../Data/lightminoil.csv')
    #df_lmol.plot(x='Time', y='ENV')

```

```

#!/usr/bin/env python3
# -*- coding: utf-8 -*-
"""
Created on Tue Mar 04 13:28:27 2025

@author: Aditya K. Rao
@github: @adirao-projects
"""

import pandas as pd

def load_data(name):
    df = pd.read_csv(f'{name}.csv', names=['Time', 'ENV', 'Q', 'I'], skiprows=2)
    return df

def get_uncert(df):
    df_ = df.where(df['Time']<0)
    df['uENV'] = df_[['ENV']].max()-df_[['ENV']].min()
    df['uQ'] = df_[['Q']].var()
    df['uI'] = df_[['I']].var()

    return df

def plot(df, name):
    sig = df.plot(x='Time', y='ENV', legend=False, grid='on',
                  title='Raw Pulse Data', color=['black'],
                  xlabel='Time (s)', ylabel='Observed Voltage (V)',
                  figsize=(14, 8)).get_figure()

    dif = df.plot(x='Time', y='I', legend=False, grid='on',
                  title='Raw Pulse Data', color=['black'],
                  xlabel='Time (s)', ylabel='Difference Voltage (V)',
                  figsize=(14, 8)).get_figure()

    sig.savefig(f'{name}sig')
    dif.savefig(f'{name}dif')

if __name__ == '__main__':
    # LMO : 56, 57, 58
    # HMO : 60, 61, 62
    # H2O : 64, 65, 66
    scope_files = ['56', '57', '58', '60', '61', '63', '64', '65', '67', '68']
    #scope_files = ['46', '47', '48', '49']
    for scp in scope_files:
        df = load_data(f"scope_{scp}")
        plot(df, scp)

```

